



Adding Monitors

Additional physiological parameter monitors can be added at any time throughout the duration of the Anaesthesia Record. The newly added monitor will only be documented from that time forward.


1. Click on **Document** in the top menu bar of Anaesthesia.
2. Click **Monitors**.



3. The Select Monitors window opens.



4. Click on the parameter you want to add.
5. Click **OK**.
6. The new parameter is added to the gases/ monitors list. If the gas or monitor is not immediately visible, scroll up or down.



Note: Parameters/ monitors chosen must match the 'label' the physiological monitoring device has attached to it. For example: ABP compared to ART; Temp compared to Temp naso; Sats compared to Sats po.